

OMLA-B

Угловые кронштейны



Преимущества серии

- Высокопрочный материал из алюминиевого сплава 6061-T4 с точной обработкой и черным анодированием.
- Минимизирует воздействие на оптическую систему путем уменьшения вибраций поверхности.
- Обеспечивает соединение между стандартными ручными позиционерами OMTOOLS для создания XZ и XYZ конфигураций.






Описание серии

Угловые кронштейны OMTOOLS подходят для ручных позиционеров и часто используются для их ортогональной установки. Линейный XY позиционер и линейный позиционер могут быть объединены в XYZ конструкцию с угловым кронштейном. Или два линейных позиционера могут быть собраны в XZ конфигурацию.

Отличающиеся параметры моделей





Модель	Размеры, мм	Совместимость
OMLA01B	65 x 65 x 90	SPYB25, SPYB50, SPYD65
OMLA02B	65 x 65 x 120	SPYB75, SPYB100
OMLA03B	40 x 40 x 58	SPYD40
OMLA04B	60 x 60 x 80	SPYD60
OMLA05B	80 x 80 x 105	SPYD80
OMLA06B	70 x 100 x 125	SPYD100
OMLA07B	70 x 125 x 150	SPYD125

Чертеж OMLA01B








DRAWING PROJECTION		 	OMTOOLS	
REV	A0	OMLA01B		
DRAWN	TODD			
DATE	2021/5/18	WEIGHT	SCALE	
FOR INFORMATION ONLY NOT FOR MANUFACTURING PURPOSES		0.23kg	1:2	

Чертеж OMLA02B








DRAWING PROJECTION					
REV	A0	OMLA02B			
DRAWN	TODD				
DATE	2021/5/20	WEIGHT	SCALE		
FOR INFORMATION ONLY NOT FOR MANUFACTURING PURPOSES		0.28kg	1:2		

Чертеж OMLA03B






DRAWING PROJECTION		OMTOOLS	
REV	A0		
DRAWN	TODD	OMLA03B	
DATE	2021/5/20	WEIGHT	SCALE
FOR INFORMATION ONLY NOT FOR MANUFACTURING PURPOSES		0.07kg	1:2

Чертеж OMLA04B







DRAWING PROJECTION			
REV	A0		
DRAWN	TODD	OMLA04B	
DATE	2021/5/20	WEIGHT	SCALE
FOR INFORMATION ONLY NOT FOR MANUFACTURING PURPOSES		0.19kg	1:2

Чертеж OMLA05B







DRAWING PROJECTION					
REV	A0	OMLA05B			
DRAWN	TODD				
DATE	2021/5/20	WEIGHT	SCALE		
FOR INFORMATION ONLY NOT FOR MANUFACTURING PURPOSES		0.42kg	1:2		

Чертеж OMLA06B



DRAWING PROJECTION			
REV	A0		
DRAWN	TODD	OMLA06B	
DATE	2021/5/20	WEIGHT	SCALE
FOR INFORMATION ONLY NOT FOR MANUFACTURING PURPOSES		0.54kg	1:2

Чертеж OMLA07B



DRAWING PROJECTION					
REV	A0	OMLA07B			
DRAWN	TODD				
DATE	2021/5/20	WEIGHT	SCALE		
FOR INFORMATION ONLY NOT FOR MANUFACTURING PURPOSES		0.78kg	1:2		