







φ2.01mm [0.079in]
2.5mm [0.1in] DEEP
ALIGNMENT PIN FEATURES
2 PLACES


NYLON-TIPPED SET SCREW FOR
SECURING OF MOUNTED OPTIC
2.0mm [5/64in] HEX



MOUNTING DOUBLE-BORE ACCEPTS 25.4mm [1.00in]
MINIMUM 3.0mm[0.12in] THICK OPTICS
MAXIMUM 7mm[0.28in] THICK OPTICS
6.3mm[0.25in] DEEP
22.8mm [0.9in]CLEAR APERTURE

M6X0.25 ADJUSTER SCREWS
PROVIDES +/-4° TOTAL ANGULAR DISPLACEMENT
3 PLACES
2.0mm [5/64in] HEX



DRAWING PROJECTION 		OMTOOLS	
REV	A0	OMBT25AN 不锈钢光学调整架	
DRAWN	TODD		
DATE	2022/9/2	WEIGHT	SCALE
FOR INFORMATION ONLY NOT FOR MANUFACTURING PURPOSES		0.24kg	1:1