



M4 x 0.7 TAPPED HOLE
6.4mm [0.25in] DEEP


SELF-CENTERING LENS CELL
MOUNTS AND POSITIONS OPTICS
BETWEEN $\varnothing 3.8\text{mm}$ [$\varnothing 0.15\text{in}$] AND $\varnothing 45.0\text{mm}$ [$\varnothing 1.77\text{in}$]
ABOUT CAGE SYSTEM'S OPTICAL AXIS




M3 LOCKING SET SCREW
1.5mm [0.08in] HEX
8 PLACES



NON-MARRING
OPTIC CLAMPING FINGER
3 PLACES



DRAWING
PROJECTION



OMTOOLS

REV

A0

OMDH45BN2

DRAWN

TODD

DATE

2023/3/10

WEIGHT

SCALE

FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES

0.15kg

1:1.2