



SM2 (2.035in-40) SERIES INTERNAL THREAD  
ACCEPTS 50.8mm [ 2.00in] MAXIMUM 7.5mm[0.3in] THICK OPTICS  
2 GYKH50A RETAINING RINGS INCLUDED  
2.5mm [0.1in] THICK  
48.3mm [ 1.90in] CLEAR APERTURE




Ø 6.0mm [ 0.24in] THRU HOLE  
FOR USE WITH CAGE RODS  
4 PLACES



LOCKING FOR 6mm ( 0.24in ) RODS  
1.5mm(0.06in) HEX  
8 PLACES



|                                                                                                          |           |                |       |
|----------------------------------------------------------------------------------------------------------|-----------|----------------|-------|
| DRAWING PROJECTION  |           | <b>OMTOOLS</b> |       |
| REV                                                                                                      | A0        | LBB50B         |       |
| DRAWN                                                                                                    | TODD      |                |       |
| DATE                                                                                                     | 2020/12/1 | WEIGHT         | SCALE |
| FOR INFORMATION ONLY<br>NOT FOR MANUFACTURING PURPOSES                                                   |           | 0.05kg         | 5:6   |