






SM1 (1.035in-40) SERIES INTERNAL THREAD
 ACCEPTS 25.4mm [1.00in] MAXIMUM 6.9mm[0.27in] THICK OPTICS
 1 GYKH25A RETAINING RINGS INCLUDED
 2.1mm [0.08in] THICK
 23.0mm [0.90in] CLEAR APERTURE



DRAWING PROJECTION			OMTOOLS	
REV	A0	OMXY30B		
DRAWN	TODD			
DATE	2023/7/19	WEIGHT	SCALE	
FOR INFORMATION ONLY NOT FOR MANUFACTURING PURPOSES		0.01kg	2:1	