







Ø 6.0mm [0.24in] THRU HOLE
FOR USE WITH CAGE RODS
2 PLACES

SM1 (1.035in-40) SERIES INTERNAL THREAD
ACCEPTS 25.4mm [1.00in]
MAXIMUM 6.1mm[0.24in] THICK OPTICS
1 GYKH25A RETAINING RINGS INCLUDED
2.0mm [0.08in] THICK
23.0mm [0.90in] CLEAR APERTURE



DRAWING PROJECTION 		OMTOOLS	
REV	A0	LBB17	
DRAWN	TODD		
DATE	2020/7/24	WEIGHT	SCALE
FOR INFORMATION ONLY NOT FOR MANUFACTURING PURPOSES		0.01kg	3:2