



| 型号        | A/mm | B/mm | C/mm | D/mm | E/mm | F/mm | G/mm |
|-----------|------|------|------|------|------|------|------|
| HGZT0606D | 600  | 600  | 800  | 100  | 380  | 380  | 120  |
| HGZT0906D | 900  | 900  | 800  | 100  | 680  | 380  | 120  |
| HGZT1008D | 1000 | 800  | 800  | 100  | 780  | 580  | 120  |
| HGZT1208D | 1200 | 800  | 800  | 100  | 980  | 580  | 120  |
| HGZT1509D | 1500 | 900  | 800  | 150  | 1200 | 600  | 200  |
| HGZT1510D | 1500 | 1000 | 800  | 150  | 1200 | 700  | 200  |
| HGZT1512D | 1500 | 1200 | 800  | 150  | 1200 | 900  | 200  |
| HGZT1609D | 1600 | 900  | 800  | 150  | 1300 | 600  | 200  |
| HGZT1812D | 1800 | 1200 | 800  | 200  | 1500 | 900  | 200  |
| HGZT1815D | 1800 | 1500 | 800  | 200  | 1500 | 1200 | 200  |
| HGZT2010D | 2000 | 1000 | 800  | 200  | 1720 | 700  | 200  |
| HGZT2012D | 2000 | 1200 | 800  | 200  | 1720 | 900  | 200  |
| HGZT2412D | 2400 | 1200 | 800  | 200  | 2120 | 900  | 200  |
| HGZT2415D | 2400 | 1500 | 800  | 200  | 2120 | 1200 | 200  |
| HGZT3010D | 3000 | 1000 | 800  | 200  | 2740 | 700  | 200  |
| HGZT3012D | 3000 | 1200 | 800  | 200  | 2740 | 900  | 200  |



|                                                                                                          |           |                  |       |
|----------------------------------------------------------------------------------------------------------|-----------|------------------|-------|
| DRAWING PROJECTION  |           | <b>OMTOOLS</b>   |       |
| REV                                                                                                      | A0        | HGZT系列阻尼平台 (4支腿) |       |
| DRAWN                                                                                                    | TODD      |                  |       |
| DATE                                                                                                     | 2022/3/24 | WEIGHT           | SCALE |
| FOR INFORMATION ONLY<br>NOT FOR MANUFACTURING PURPOSES                                                   |           | kg               | 1:20  |