



SM1 (1.035"-40)
EXTERNAL THREADS

Ø 16.1mm [0.634in]
CLEARANCE HOLE

DOUBLE BORE
PROVIDING 3 POINTS CONTACT
WITH NYLON-TIPPED SETSCREWS



M3x0.5 NYLON-TIPPED
SECURING SETSCREW
1.5mm[0.06in] HEX

DRAWING PROJECTION 		OMTOOLS	
REV	A0	FZ1SM16	
DRAWN	TODD		
DATE	2020/7/23	WEIGHT	SCALE
FOR INFORMATION ONLY NOT FOR MANUFACTURING PURPOSES		0.01kg	2:1